

- 1) Turn on power supply to 150V
- 2) Adjust head so electron path is circular & not helical
- 3) Measure left & right side of beam, calculate center & radius.

The purpose of this lab is to confirm the established value of e/m for electrons through experimentation with Helmholtz coils.


D1

Measurements:

4.5 to 5.5 cm

$r_{\text{center}} = 10.0 \pm 0.5 \text{ cm}$

$r = 5.5 \pm 0.5 \text{ cm}$



D2

- 1) Max A: 2.0 A
Min A: 0.7 A

(A) ± 0.1	(cm) ± 0.2	(cm) ± 0.2	(cm) ± 0.2
Amperage	left	right	diameter
0.7	3.25	15.75	12.5
1.0	5.25	14.75	9
1.35	5.75	13.4	7.65
1.7	6.4	12.6	6.2
2.0	6.8	12.4	5.6

- 1) Increase & decrease voltage, repeat measuring the beam across 5 points including min & max

D3

0.7 A

(V) $\pm 1\%$	(cm) ± 0.2	(cm) ± 0.2	(cm) ± 0.2
Voltage	left	right	diameter
97	4.7	13.6	8.9
105	4.4	14.3	9.9
115	4.0	14.7	10.7
125	3.6	15.1	11.5
135	3.4	15.3	11.9
145	3.1	15.8	12.7
150	2.9	15.7	13.0

$$\frac{1}{2}mv^2 = eV$$

$$\vec{F} = -e(\vec{v} \times \vec{B})$$

$$evB = mv^2/R$$

$$\frac{e}{m} = \frac{2V}{(BR)^2}$$

$$B = \frac{8\mu_0 NI_c}{5.15}$$

↓

$$\frac{e}{m} = \frac{2V}{(8\mu_0 NI_c / 5.15)^2 R^2}$$

m = mass of electron

v = speed of electron

\vec{v} = velocity

\vec{B} = magnetic field

e = magnitude of charge of electron

$$\frac{e}{m} \approx 1.75882017 \pm 0.00000007 \times 10^{11} \frac{C}{kg}$$

$$\mu_0 = 4\pi \times 10^{-7} T \cdot m/A$$

I_c = current

N = turns in coil = 130

$r = 0.158 \pm 0.005 m$

R = radius of beam

B = mag field

V = m. voltage

1) derive the relations between $R, I_c, V, e/m$ to ultimately calculate e/m

Varying current:

$$\frac{1}{R} = \frac{8\mu_0 N}{5.15} \sqrt{\frac{e/m}{10V}} I_c + C$$

Varying Voltage:

$$R = \frac{5.15}{8\mu_0 NI_c} \sqrt{\frac{10}{e/m}} \sqrt{V} + C$$

$$\delta \frac{1}{R} = \frac{1}{R^2} \delta R$$

$$\delta \sqrt{V} = \frac{1}{2\sqrt{V}} \delta V$$

- 1) Connect calibration to CH1 of the DSO
- 2) measure Volts/DIV, time/DIV
- 3) estimate error in DSO
- 4) Measure period & peak-to-peak of the calibration wave.


The purpose of this lab is to familiarize myself with the new generator and the DSO, and also to explore how best to work.

$$V/DIV = \cancel{500mV} \quad 500mV \quad ms/DIV = 250$$

$$error = \frac{1}{5} DIV$$

$$period = \frac{16}{4} DIV = 1000ms \pm 10ms$$

$$peak-to-peak = 6 DIV = 3V \pm 100mV$$



$$\frac{V}{DIV} = 2.00V$$

$$\frac{ms}{DIV} = 5.00ms$$

$$peak-to-peak = 6 DIV = 12.0 \pm 0.4 V$$


$$period = \frac{16.6 DIV}{5} = \cancel{16.6ms} \quad 16.6 \pm 0.2 ms$$

$$DMM \text{ reads } 4.27 \pm 1\%$$

$$V_{pp} = 2\sqrt{2} V_{rms}$$

$$V_{pp} = 12.077, \text{ this is consistent with our oscilloscope reading.}$$


60 Hz:




1) add wave generator to 412 input, these patterns are XY plots of sine wave vs. wave generator.

120 Hz:

$$V = 120.043 V$$



90 Hz: The ratios of the frequencies need to be a simple ratio.



$$\frac{314 \text{ ns} \pm 5}{30} = 10.47 \pm 0.2$$

$$10.5 \pm 0.2 \text{ ns}$$

1) measure waveform of sound, calculate frequency of it.